





TRI-Vent<sup>o</sup>

**TRI**  
surgical kit

**Through Research Innovative**

# TRI-Vent<sup>®</sup> Surgical Kit



- |                                             |                                                              |
|---------------------------------------------|--------------------------------------------------------------|
| <b>1</b> Extra                              | <b>16</b> PP – Paralleling Tool                              |
| <b>2</b> Extra                              | <b>17</b> PP – Paralleling Tool                              |
| <b>3</b> DE – Drill Extender                | <b>18</b> TID 2,5 – Implant Driver-Short, 2,5 mmD Hex        |
| <b>4</b> Pilot Drill – Short, 2,3 mmD       | <b>19</b> TIDL 2,5 – Implant Driver-Long, 2,5 mmD Hex        |
| <b>5</b> Pilot Drill – Long, 2,3 mmD        | <b>20</b> TIDD 2,5 – Implant Driver, 2,5 mmD Hex – Handpiece |
| <b>6</b> Step Drill – Long, 2,5/2,8 mmD     | <b>21</b> THD 1,25 – Hex Driver-Short, 1,25 mmD              |
| <b>7</b> Step Drill – Short, 2,5/2,8 mmD    | <b>22</b> THDL 1,25 – Hex Driver-Long, 1,25 mmD              |
| <b>8</b> Step Drill – 2,8/3,5 mmD           | <b>23</b> THDD 1,25 – Hex Driver, 1,25 mmD – Handpiece       |
| <b>9</b> TAP – Cortical Bone Drill 3,75 mmD | <b>24</b> TART – Removal Tool for Abutments                  |
| <b>10</b> Step Drill – 3,0/3,5 mmD          | <b>25</b> Extra                                              |
| <b>11</b> Step Drill – 3,3/3,9 mmD          | <b>26</b> Extra                                              |
| <b>12</b> TAP – Cortical Bone Drill 4,1 mmD | <b>27</b> Extra                                              |
| <b>13</b> Step Drill – 3,5/3,8 mmD          | <b>28</b> Ratchet                                            |
| <b>14</b> Step Drill 3,8/4,5 mmD            | <b>29</b> Surgical Driver                                    |
| <b>15</b> TAP – Cortical Bone Drill 4,7 mmD |                                                              |

